

ENTity SDXL Flexible Fiberoptic Laryngoscope


Designed to empower ENT clinicians with the visibility & mobility they need



Slimline 2.7mm shaft makes the ENTity SDXL the preferred choice for narrow upper airway passages

The ENTity SDXL Flexible Fiberoptic Laryngoscope frees you to perform nasal endoscopy and laryngoscopy virtually anywhere. Our next generation patented, integrated LED light source eliminates cumbersome external light sources while providing exceptional visibility.

The ENTity SDXL now features a *maximum brightness* 18% *higher than the previous generation*. Additionally, the new generation LED system in the endoscopes provide a vibrant white light targeted in the 5500K color temperature range for superior tissue color differentiation.



Optilux^{*}

Optim's unique LED light source to fiber light guide coupling technology, Optilux™ in the ENTity SDXL brings unmatched illumination to the end of the scope, where it counts.

18% Increase in maximum brightness

Compact, portable, all-in-one design

- No external light source or lamp replacement
- Hybrid aluminum handle
- Fully sealed and completely immersible
- Bring the scope to the patient
- Use in the office / bedside / hospital

Better economics

- See more patients per day, increase revenue
- Reduced capital expenses and ongoing maintenance
- No lamp or light guide replacements

Superior lighting

- Patented LED light coupling
- Enables clear visibility of the true vocal cords
- Variable light intensity

Specifications

Flexible Length	30cm (11.8")
Overall Length	53cm (20.9")
Shaft Diameter	2.7mm
Weight	11.6 oz
Articulation Range	110 degrees up / down
Depth of Field	5mm to 50mm
Field of View	85 degrees
Battery	Lithium Ion



64 Technology Park Road, Sturbridge, MA 01566 Toll-free: 800.225.7486 Worldwide: 1.508.347.5100

Fax: 508.347.2380

Email: sales@optim-llc.com Website: www.optim-llc.com

