

The Sonesta line of chairs

Top quality | Versatile | User friendly | Infinite positioning



Sonesta manufactures and sells tables and chairs for urology, gynecology, urodynamic diagnostics and upper GI procedures

Next Generation Sonesta

SONESTA S2



Specifications

Length with extension cushions	193 cm	76"
Length without	143 cm	56"
Width	77 cm	30"
Height with seat cushion (min)	40 cm	16"
Height with seat cushion (max)	105 cm	41"
Lift capacity	210 kg	465 lbs
Wheel diameter (Locking)	10 cm	4"
Net weight	100 kg	220 lbs

Electrical

120/230v 24v DC Output
Meets EMC (Electromagnetic Compatibility) for IEC 60601
All Products carry the CE mark.

Features:

- Versatility
- Precise positioning
- Two motors are offered for height, and back adjustments
- Easy access for the elderly and handicapped
- Engineered Ergonomics
- Infinite adjustments as low as 16" (40cm) to 41" (105cm) high
- Light weight and easy to move
- Supporting 465 weight capacity
- UL Certified

Standard Accessories:

2 Arm rests	525-6339-2
2 Leg rests	525-6336
2 Leg rest holders	525-6338
1 Paper towel holder	525-6326
1 Two motor foot control	524-6353
1 Funnel	524-6325-5

Next Generation Sonesta

SONESTA S3



Specifications

Length with extension cushions	193 cm	76"
Length without	143 cm	56"
Width	77 cm	30"
Height with seat cushion (min)	40 cm	16"
Height with seat cushion (max)	100 cm	39"
Lift capacity	210 kg	465 lbs
Wheel diameter (Locking)	10 cm	4"
Net weight	125 kg	276 lbs

Electrical

120/230v 24v DC Output
Meets EMC (Electromagnetic Compatibility) for IEC 60601
All Products carry the CE mark.

Features:

- Locking Mechanism for one touch movement for precise positioning
- Versatility
- Three motors are offered for height, seat and back adjustments
- Easy access for the elderly and handicapped
- Engineered Ergonomics
- Infinite adjustments as low as 16" (40cm) to 39" (100cm) high
- Light weight and easy to move
- Supporting 465 weight capacity
- UL Certified

Standard Accessories:

2 Arm rests	525-6339-2
2 Leg rests	525-6336
2 Leg rest holders	525-6338
2 Fold up footrests	525-6331
1 Seat cushion without cut out	524-6365
1 Cushion with cut out	525-6366-3
1 Paper towel holder	525-6326
1 Three motor foot control	524-6357
1 Funnel	524-6325-5
1 Neck pillow	525-6346



About Next Generation Sonesta:

Sonesta is a global company focusing on quality, function and comfort. We have a large install base in the world and have now created an improved chair.

Our goal was to develop the next generation Sonesta chair, while at the same time maintain the "Sonesta" branded look and to enhance the feel with improved human interaction. We wanted to give the patient, the doctors and nurses a more holistic pleasurable user experience.

We are proud to introduce the Next Generation Sonesta chair. The Next Generation Sonesta chair is the most versatile chair on the market and the ultimate patient positioning chair with infinite positioning.



Procedures commonly used in the Sonesta chair:

- Urodynamic
- Pelvic Muscle Rehabilitation
- Anorectal Manometry
- Pelvic Exams
- Cystoscopy
- Ultrasounds
- Vasectomy
- Colposcopy
- Biopsy's
- Colonoscopy
- Gynecology procedures



Features & Benefits of the Next Generation Sonesta chair:

Weight Capacity - 465 pounds.

New soft ergonomic and low profile cushions with improved foam - To enhance patient comfort and feel.

Fold up arm rests - Easy access to the chair from each side as well as a safety protection when chair is in a table position.

Height adjustable head rest - For patient comfort, and easily adjustable.

Vacuum formed back cover - Protecting and covers the frame and other structural parts, for a calmer and less mechanical feeling. Easier to clean.

Neck Pillow - Easier to adjust height using magnets for a more premium feel.

New handlebar - To enhance moving of the chair, with an intergraded paper towel holder.

Paper towel holder - User friendly easier to change paper roll.

Footrests - Has improved with a lever, to ease unlocking to insure safety for end users. Plastic covered magnet added to keep the foot rest in place when folded up.

New ergonomic leg rest - With improved comfort and asymmetric shapes makes it possible to reverse the leg rest on usage purpose.

Fluoroscopy Procedure Table

SONESTA 6210



Specifications

<i>Lengths without foot rest/leg extension</i>	156 cm	63"
<i>Length with foot rest/leg extension</i>	196 cm	78"
<i>Width (Total)</i>	99 cm	39"
<i>Height with seat cushion (min)</i>	50 cm	20"
<i>Height with seat cushion (max)</i>	130 cm	51"
<i>Imaging area without foot rest/leg extension (length)</i>	156 cm	63"
<i>Imaging area with foot rest/leg extension (length)</i>	196 cm	78"
<i>Weight capacity</i>	250 kg	551 lbs

Electrical

120/230v 24v DC Output
 Meets EMC (Electromagnetic Compatibility) for IEC 60601
 All Products carry the CE mark.

A versatile examination table with C-arm capabilities

The 6210 fluoroscopy procedure table offers the physician standing, seated, and supine positioning capabilities. This versatile table easily adjusts from/to standing, seated, and supine positions via a hand held control, thus creating a table that is easily accessible to patients of any age. The side-mounted cantilevered design offers the physician anterior/posterior, lateral and oblique views without obstruction. For added convenience the table is equipped with a battery backup for short-term use in case of power failures or in emergency situations. The superb c-arm access and a wide range of positioning possibilities, make the 6210 the ideal table for any fluoroscopy procedure.

Positioned to be the best:

- Three motors offer unparalleled positioning/access options
- Radiolucent carbon fiber seat and back for head-to-toe imaging
- 4 Double wheel rolling casters offer unencumbered maneuverability



Accessories:

1 Hand control	524-6356
1 Plastic funnel w/drainage	524-6325-9
1 Arm rest holder	524-6338-4
1 Arm rest	524-6021-6
2 Leg rest holder	524-6338
2 Heavy duty Rissler Leg rest	524-6336-4
2 Fold up foot rest	524-6331
1 Side Rail	524-6329

An ADDvise Group Company



SONESTA MEDICAL AB
 TEGELUDDSVÄGEN 76
 115 28 STOCKHOLM, SWEDEN
 TEL: +46 8 502 572 80
 INFO@SONESTAMEDICAL.SE
 WWW.SONESTAMEDICAL.SE

SONESTA MEDICAL INC.
 2 DEBUSH AVENUE UNIT C3
 MIDDLETON, MA 01949, USA
 TEL: +1 630-519-3450
 INFO@SONESTAMEDICAL.SE
 WWW.SONESTAMEDICAL.SE